

Jégu T¹, Thion C¹, Murphy A², Monjaret F¹, Michaelis J³, Plesshoff S³, Mykhailiuk K, Weigel W³, Berthuy O¹, Cantlon J² and Tourniaire G¹


¹Cellenion SASU, Lyon, France, ²SCIENION US Inc., Tempe, AZ, United States, ³SCIENION AG, Berlin, Germany

www.cellenion.com – info@cellenion.com

INTRODUCTION


Single-cell RNA sequencing (scRNA-seq) methods must combine sensitivity of capturing transcript diversity within a cell with an accuracy that reflects the relative proportions of the transcripts in a cost and time-efficient manner. This preliminary work shows that the cellenONE[®] platform, providing both single cell isolation and automated nanoliter dispensing, drastically reduces library prep cost and handling time by miniaturizing an existing scRNA-seq commercial kit, from 5- μ l \rightarrow 100-nl working volume.

WORKFLOW



METHODS

Nanoliter liquid dispensing & single-cell isolation with **cellenONE[®] X1**



Library preparation reagents from **QIAseq UPX 3' Transcriptome kit**


Nanoliter volume reactions in **cellenCHIPS[®]** (4x96-wells in a microscope slide footprint)

+ Human embryonic kidney cells = **HEK**

or Chinese hamster ovary cells = **CHO**

RESULTS: Sensitivity

Thousands of transcripts captured and genes detected per cell:



With **100 nl RT reaction**:

- ✓ Up to **23,000 captured mRNA molecules** and **2,800 detected genes** per unique cell (average 8,000 and 1,200 respectively)
- ✓ **116,000 captured transcripts** and **6,500 detected genes** in "Multiple" (<10) cells
- ✓ **Comparable to published QIAGEN results in original microliter-volume protocol** and to other scRNA-seq technologies (Ziegenhain, 2016, Mol Cell DOI:10.1016/j.molcel.2017.01.023!)

RESULTS: Contamination

A clean signal with low background:




- ✓ **Low to null** cross-sample contamination
- ✓ **Low to null** library background noise

RESULTS: Robustness

Reproducibility highlights transcriptome differences:

- ✓ **Replicate** cellenCHIPS (■ and ▲) do not discriminate from one another
- ✓ Exp. 1 cell transcriptomes (●) significantly discriminate from Exp. 2 transcriptomes (■ and ▲)



- ✓ Hamster cells (■, ▲ and ●) always neatly discriminate from human cells (■, ▲ and ● resp.)

CONCLUSION

- ✓ Miniaturized Qiagen's UPX 3' Transcriptome kit for use in cellenONE[®] with cellenCHIP[®] and obtained **scRNA-seq libraries for less than 1.75 € / 2 USD per cell**
- ✓ Data quantitatively aligns with other scRNA-seq technologies
- ✓ Current work focuses on improving robustness / sensitivity:

Watch out for an Application Note and full protocol details, to be released soon!

Working solutions (scaling quantities to 30-60 and 50-100 nl reaction volumes for cell lysis and RT reaction mixes, respectively) were prepared, and incubations were performed, according to manufacturer's instructions and data was processed using GeneGlobe, a dedicated Qiagen pipeline.