



cellenONE® Sterile Condition Operation

Cellenion SASU, 60 Avenue Rockefeller, BioSerra2, 69008 Lyon, France

Abstract

Single cell methods require extreme cleanliness to avoid exogenous contamination, as a minute microbial contamination can ruin a cell cloning procedure and jeopardize a single cell sequencing approach. We evaluated axenic (germ free) conditions with cellenONE® comparing two systems: a regular cellenONE® X1 and a cellenONE® X1 BSC (version of the instrument mounted in a BSC environment).




Figure 1. Sterilisation test workflow

Materials and methods


Sterilization experiments consisted of 3 conditions (Fig. 1A), where the dispense was performed for each phase (Fig. 1B); 1) Before contamination: in **routinely clean** (H₂O, EtOH, detergent) cellenONE® systems. 2) Contamination: in **intentionally contaminated** systems. Environment generated by aspirating and dispensing a concentrated suspension of viable *E. coli* bacteria and *A. niger* fungal spores with the cellenONE® nozzle. 3) After decontamination: in **sterilized** systems.

1. Drops were **cellenONE® dispensed in sterile growth media** (LB, DMEM and PDB for bacterial, human and fungal cells, respectively, 100 nL in 250 µL) and **PCR water** (2 µL in 38 µL).
2. Inoculated media were incubated (48h, 37°C) then spotted on **agar plates**, to assess for **culturable** bacterial and fungal load.
3. Inoculated PCR water was used as template for **real-time PCR** (qPCR) targeting universal bacterial 16S ribosomal genes and fungal ITS genes, to assess **total** microbial load.



Results and discussion

In **routinely clean** cellenONE® X1, 1 out of 48 wells had a bacterial contamination, none in X1 BSC.



All **positive control** wells (intentionally heavy microbial contamination) were contaminated as expected and flushing the nozzle with water was not sufficient to remove the contamination.


cellenONE® X1 and cellenONE® X1 BSC **became germ-free** again after sciCLEAN8 and Sterilization tasks, respectively.



In the X1, the sterilization task* induces temporary sterile conditions. The X1 BSC is an operational germ-free environment such that detergent-based washing, sciCLEAN8, is sufficient to prevent contamination.



B: bacteria, **F:** fungi,
Cult. and **Tot.:** culturable and total
microbial load
+++ : positive control
++ : major contamination (>1/12 wells)
+ : minor contamination (1/12 wells)
- : no contamination



			Cult.	Tot.	Cult.	Tot.
Before		B	+	-	-	-
		F	-	-	-	-
Contamination	Microbial suspension	B	+++	+++	+++	+++
		F	+++	+++	+++	+++
	After flushing nozzle	B	++	++	+	-
		F	++	-	-	+
After decontamination	sciCLEAN8	B	+	-	-	-
		F	+	-	-	-
	Sterilization Task*	B	-	-	-	-
		F	-	-	-	-

** A sequence of nozzle washing with 0.5 % sodium hypochlorite, 3% hydrogen peroxide and 70% ethanol*

Contact Us

Cellenion SASU
60 Avenue Rockefeller
Bioserra 2, 69008
LYON France
Tel: +33 986 48 70 70
contact@cellenion.com
www.cellenion.com

SCIENION GmbH
Volmerstr. 7b
D-12489 Berlin
Germany
Tel: +49 (0)30 6392 1700
support@scienion.com
www.scienion.com

SCIENION US, Inc
4405 E. Baseline Road Suite
#123
Phoenix, AZ. 85042
United States
Tel: +1 (888) 988-3842
USsalesupport@scienion.com

SCIENION (UK) Ltd
2000, Lakeside North Harbour
Western Road, Portsmouth
PO6 3EN
United Kingdom
Tel: +44 (0) 1243 88 71 65
support@scienion.com