

Single cell isolation with
cellenONE[®]
for sequencing applications

with integrated nanoliter dispenser



scRNA-Seq




scWGS



scATAC-Seq



scMethyl-Seq



scMS



scYour-Seq


LOW VOLUME, HIGH PRECISION & RECOVERY SINGLE CELL ISOLATION
FOR SEQUENCING... AND MUCH MORE

About the technology


cellenONE, an automated single cell dispensing system based on patented piezo-acoustic technology, allows precise cell deposition on a wide range of microplates (96, 384, 1536) and microwell substrates.

Most dispensing and microfluidic technologies follow Poisson distribution, which leads to multiple cells per position, low efficiency and biased data. cellenONE uses software-integrated visual feedback to ensure only single cells are deposited in every position.


Each drop generated with cellenONE can contain:



(A) no cell



(B) one cell,





(C) multiple cells.

cellenONE only dispenses single-cell containing drops (B) directly into microplates or microwell chips of your choice.


All remaining drops are dispensed into a recovery tube, resulting in none of the cells being lost and leaving the possibility to reprocess those non-isolated cells later.

Samples are analyzed live and pictures recorded for post processing analysis



Single cell deposited in 50 μm microwells

Thanks to visual feedback, only single cells are dispensed








NOW also available with fluorescence detection module

cellenONE can now be ordered with optional fluorescence detection setup which can be used for:

- Isolation of live single cells (Live/Dead)
 - Remove dead cells from scRNA-seq libraries*
- Isolation of very rare single cells
 - Only isolate single cells of interests*
- Cloning best expressers
 - Boost clonal recovery and expressers levels*

Optional fluorescence detection comes in two levels:

One channel	FITC
Four channels	DAPI
	FITC
	PE
	Cy5



Four arrays of 5x5 single HEK293T cells labeled with DAPI, CellTracker Green, Orange and Deep Red sorted using 4-channels setup.

A mixture of four cell populations were labeled with Hoechst 33342 (DAPI channel), CellTracker Green (FITC channel), CellTracker Orange (Cy3 channel) and CellTracker Deep Red (Cy5 channel).


25 single cells from each cell population were successfully isolated on microscope slide with the fluorescence detection module of cellenONE.

Sorting Features

Live / Dead


Very rare cells

Only isolate single cells of interest



Up to 100% single cells

Outstanding accuracy



Results from 5x100 positions filled with single cells from four different cells samples. Single cells rate is indicated in pink, up to 99% for A549 cells; in white are empty positions and in grey, multiple cells positions.


Unlike other technologies (microdroplet, microfluidic traps), doublets can be completely excluded and only single cells are dispensed.

Open and versatile platform

Thanks to high definition optics, cells are individually identified and isolated. These are selected according to their size, geometry and now potentially using fluorescent markers (see page 3). Once selected, cells are dispensed into microplates such as 96, 384, 1536wp or into microwells such as ICELL8™ chip.

Successfully isolated cells and particles so far:

- **cell lines** such as CHO, hybridoma, HEK293T, HeLa, A549, PC3, H1975, HepaRG, Jurkat
- **primary cells** such as PBMC (including B- and T- cell fractions), fibroblasts, keratinocytes, melanocytes, cardiomyocytes, HUVEC, neural stem cells
- **nuclei** from cell line, fresh frozen (FF) and formalin-fixed paraffin-embedded (FFPE) tissue slices
- **microbeads** like PMMA microbeads from 2 to 30µm diameters, TOYOPEARL®, PS microbeads



Single cells from dissociated lung cancer spheroids successfully isolated onto a microscope slide. Every position contains one single cell.


Key Features

cellenONE®

No doublets

Unseen recovery

Suitable for most cell types



Cells and particles from 2 to 70 µm

System Set up Workflow



Cell Sample Preparation

Place sample in holder
Aspirate sample

- Prepare a minimum volume of 2 μL , ideally around 10 μL
- Optimal cell concentration below 200 cells/ μL




Cell Mapping

- Automated cell tracking and report generation
- Choose isolation parameters
- Select type of target (96, 384, 1536 wp. ICELL8™ chip)

System setup only takes a few minutes

Users can walk away during isolation process



A mapping report is automatically generated, it contains:


- Evaluation of ejection zone boundary
- Detected cell diameter, elongation and circularity distribution
- Average cell concentration



Automated Isolation

- Up to 2 microplates (or equivalent) per run
- 100 cells in 4 min or 1000 cells in less than 45 min*
- Walk-away and easy to operate system
- Temperature control to avoid RNA degradation

* at optimal cell concentration



Open Workflow for Library Preparation



scRNA-Seq



scWGS



scATAC-Seq



scMethyl-Seq




scYour-Seq

cellenONE is an open platform providing both single cell isolation and nanoliter dispensing.

Such versatility allows users to automate many of the crucial steps involved in the ever growing number of single cell library preparation methods available today. Moreover, ability to work in nanovolumes and microwells allows a drastic reduction in reagents consumption and associated costs.

Benefits cellenONE®

- All-in-one device
- Low reaction volumes
- Low reagent costs



Rare cells isolation from cerebrospinal fluid

Unseen recovery for low concentration samples for low volume samples




UCSF

Department of
Neurology
Weill Institute for Neurosciences

In collaboration with the Department of Neurology at Weill Institute in San Francisco, a study performed aiming to isolate single immune cells in multiple sclerosis patient's cerebrospinal fluid (CSF) samples for subsequent single cell RNA sequencing.

The main challenge associated with such sample is the very low number of cells present in these samples which can vary from only 20 to 500 cells per 20mL of patient's CSF.




To demonstrate achievable cellenONE recovery rates, a proof of concept experiment was undertaken with low concentration PBMC samples (~25 cells/ μ L).

Three samples of 30 μ L were prepared, for each sample, 4 aliquots of 5 μ L were pipetted into 4 wells of a 1536 wp and the exact number of cells was counted under microscope to calculate starting cell concentrations (average input or).

Finally, for each sample, a 5 μ L aliquot was processed using cellenONE and as many as possible single cells were isolated into a microwell chip and subsequently counted by microscopy (output or).

On average, recovery reached 95%, meaning most of the cells present in each sample were successfully isolated as single cells (number of single cells successfully isolated with cellenONE over average cell number for each sample).



Benefits
cellenONE[®]

From minute samples
(from just 2 μ L)

Recovery up to 95%

No dead volume



High throughput, versatile and low cost WGS

Scalable whole genome sequencing of 40,000 single cells identifies stochastic aneuploidies, genome replication states and clonal repertoires



Emma Laks¹ (elaks@bccrc.ca), Hans Zahn¹, Daniel Lai¹, Andrew McPherson¹, Adi Steif¹, Jazmine Brimhall¹, Justina Biele¹, Beixi Wang¹, Tehmina Masud¹, Diljot Grewal¹, Cydney Nielsen¹, Samantha Leung¹, Viktoria Bojilova¹, Maia Smith¹, Oleg Golovko¹, Steven Poon¹, Peter Eirew¹, Farhia Kabeer¹, Teresa Ruiz de Algora¹, So Ra Lee¹, M. Jafar Taghiyar¹, Curtis Huebner¹, Jessica Ngo¹, Tim Chan¹, Spencer Vattr-Watts¹, Pascale Walters¹, Nafis Abrar¹, Sophia Chan¹, Matt Wiens¹, Lauren Martin¹, R. Wilder Scott¹, Michael T. Underhill², Elizabeth Chavez¹, Christian Steidl¹, Daniel Da Costa¹, Yusanne Ma³, Robin J. N. Coope³, Richard Corbett³, Stephen Pleasance³, Richard Moore³, Andy J. Mungall³, CRUK IMAXT Consortium⁴, Marco A. Marra⁵, Carl Hansen⁵, Sohrab Shah¹ (sshah@bccrc.ca) and Samuel Aparicio¹ (saparicio@bccrc.ca)

¹BC Cancer Research Centre; ²Michael Smith Laboratories; ³Michael Smith Genome Sciences Centre; ⁴University of British Columbia

doi: <https://doi.org/10.1101/411058>


Single cell WGS study performed jointly at UBC and BCCRC aimed to analyze genome heterogeneity, mutational process and clonal evolution over tens of thousands of single cells or single nuclei from a variety of tissues. The author developed DLP+, a fully integrated pipeline from wet bench to analysis platform.

Direct Library Preparation (DLP+) is a new single cell tagmentation-based WGS method without preamplification which includes an analysis pipeline and allows clustering of thousands of cells at once



DLP protocol, adapted from Zahn et al., Nature Methods 2017

The platform was successfully used to study lineages and tissue heterogeneity according to parameters such as genome state, ploidy, copy number alteration and single nucleotide variant.




DLP+ for efficient Whole Genome Sequencing



THE UNIVERSITY OF BRITISH COLUMBIA

Direct Library Preparation Method



DLP+ protocol, adapted from Laks et al., BioRxiv 2018

Benefits

cellenONE[®]

Low volume & reagent costs

0,27 US\$/cell

High quality library

Versatile: cell and nuclei

from 5 to 80µm

Open and high throughput


1000's cell/experiments

Results

- Compatible with cell line, dissociated mouse and human tissues and nuclei from FF human breast cancer PDX tissue
- Optimized WGS protocol for nanovolume library preparation
- Open-source cloud computing pipeline for single-cell genome analytics (<http://www.cellmine.org>)
- Pseudobulk sequence from assembled single cell genomes successfully used as reference for clustering the same single cells
- Level genome ploidies, cell specific mitotic error rates, and cell specific replication status were successfully obtained from single cell data

“cellenONE technology enables the work we are doing: isolating cell and then dispensing a few nanoliters from various reagents. The system is spatially and volumetrically accurate and is able to spot each cell, one after another, into each well.”

Dr. Robin Coope,
British Columbia Cancer Genome Sciences Centre



Specifications



Technical information

Dispensing technology: piezo acoustic drop-on-demand

Dispense volume: 50-800 pL per drop

Linear drives for X/Y and spindle drive for Z

Resolution: 1 μm

Accuracy (Absolute Position): < 10 μm

Precision (Repeat Position): < 3 μm

HD camera for detection of cells or particles from 2 μm

Max. speed: 100 cm/sec

Isolation area (mm): x=180; y=120 (2 microtiter plates)

Dimensions LxWxH (mm): 650 x 700 x 1590

-> including monitor's arm L = 1300 mm

-> with door open H = 2050 mm

Weight: 205 kg

Voltage: 110 V; 220 V

Options & Software

Temperature, humidity and dew point control

2nd channel for nanoliter dispensing

Fluorescence module up to 4 channels (DAPI, FITC, Cy3 or Cy5)

Customized holders for microwells chips

Fiducial recognition and automated target alignment


Related Products & Services

cellenBEADS for calibration

cellenWASH for sterilization

cellenVIALS for recovery

cellenSERVICES for application support



cellenion SASU
60 Avenue Rockefeller
69008 LYON France
Tel: +33 986 48 70 70
contact@cellenion.com
www.cellenion.com

SCIENION AG
Volmerstr. 7b
D-12489 Berlin
Tel: +49 (0)30 6392 1700
info@scienion.com
www.scienion.com

SCIENION US, Inc.
2640 W Medtronic Way
Tempe, AZ 85281
Tel: +1 (888) 988-3842
USinfo@scienion.com
  @SCIENION_AG