Parallel measurement of transcriptomes and proteomes from same single cells using nanodroplet splitting

James M. Fulcher¹, Lye Meng Markillie¹, Hugh D. Mitchell¹, Sarah M. Williams¹, Kristin M. Engbrecht², Ronald J. Moore³, Joshua Cantlon-Bruce^{4,5}, **Johannes W. Bagnoli⁵**, Anjali Seth⁵, Ljiljana Paša-Tolić¹, Ying Zhu¹

- 1) Environmental Molecular Sciences Laboratory, Pacific Northwest National Laboratory, Richland, Washington 99354, United States
- 2) Chemical and Biological Signature Sciences, Pacific Northwest National Laboratory, Richland, WA 99352, United States
- 3) Biological Sciences Division, Pacific Northwest National Laboratory, Richland, Washington 99354, United States
- 4) Scienion AG, Volmerstraße 7, 12489 Berlin, Germany
- 5) Cellenion SASU, 60 Avenue Rockefeller, Bâtiment BioSerra2, 69008 Lyon, France

Combination of cellenONE and NanoPOTS for Nanovolume multiomics of the same single cell

- The cellenONE platform enables highly accurate single cell and nano volume dispensing
- The nanoPOTS platform provides a robust and efficient nanonvolume substrate for omics

Take Home Message:

nanoSPLITS provides multimodial expression data from the same

single cell

nanoSPLITS precisely quantifies the abundances of both mRNA transcripts and proteins
Unbiased protein detection via mass spectrometry
Very high sensitivity in both modalities
nanoSPLITS will enable greater insights into how the different modalities interact with each other

CELLENION





solutions

 The combination of both platforms forms a perfect basis for a multiomics approach (nanoSPLITS)



NanoSPLITS shows high sensitivity for proteomics and transcriptomics



Overview of the nanoSPLITS-based single cell multiomics platform.



- Lysis and single-cell disepesning using the cellenONE
- Molecule splitting by drop seperation between 2 chips

 Proteomics and transcriptomics data show similar distributions of gene ontologies within the respective datasets, e.g. for cellular components Parallel Proteomics and transcriptomics preperation from the same cell

Protein and RNA expression show poor correlation within the same single cell



0.1 0.2 0.3 0.4 0.5 0.6 0.7 0.8 0.9

C10 cells

SVEC cells

Pearson Correlation (r)

scRNAsea

scProteomics



Both modalities, scProtemoics and scTranscriptomics can easily identify cell types via their expression
 Within modality correlation are higher in proteomics data
 Across modality correlations are poor (r= 0.31-0.56)







• Combined data shows superior performance in clustering

 Identification of cell cycle derived subpopulations in C10 cells validates the high sensitivity of nanoSPLITS despite droplet splitting.

• Modalities show different top 5 markers to distuingish celltypes (C10, SVEC)

• H2-K1 can be found in both modalities as a cell type marker

Contact Us

Cellenion SASU 60 Avenue Rockefeller Bioserra 2, 69008 LYON France Tel: +33 986 48 70 70 contact@cellenion.com www.cellenion.com

SCIENION GmbHSCIENIONVolmerstr. 7b4405 E. BaselinD-12489 Berlin#123GermanyPhoenix, AZ. 8Tel: +49 (0)30 6392 1700United Statessupport@scienion.comTel: +1 (888) 98www.scienion.comUSsalessupport

SCIENION US, IncSCIENION (UK) Ltd4405 E. Baseline Road Suite2000, Lakeside North Harbour#123Western Road, PortsmouthPhoenix, AZ. 85042PO6 3ENUnited StatesUnited KingdomTel: +1 (888) 988-3842Tel: +44 (0) 1243 88 71 65USsalessupport@scienion.comsupport@scienion.com

Fulcher, J.M., Markillie, L.M., Mitchell, H.D., Williams, S.M., Engbrecht, K.M., Moore, R.J., Cantlon-Bruce, J., Bagnoli, J.W., Seth, A., Paša-Tolić, L., et al. (2022). Parallel measurement of transcriptomes and proteomes from same single cells using nanodroplet splitting.

Zhu, Y., Piehowski, P.D., Zhao, R., Chen, J., Shen, Y., Moore, R.J., Shukla, A.K., Petyuk, V.A., Campbell-Thompson, M., Mathews, C.E., et al. (2018). Nanodroplet processing platform for deep and quantitative proteome profiling of 10-100 mammalian cells. Nat. Commun. 9, 882.

©Cellenion 2022, All rights reserved