PP 13 - Enhance your 3D cellular aggregates assays' quality and reproducibility with automated single spheroids and organoids sorting and isolation from bulk populations

Steffen Cosson¹, Minjoung Jo¹, Mathieu Bennet¹, Oliver Krispin², Paul Kollhof², Martin Horn², Sébastien Clerc¹, Guilhem Tourniaire¹ ¹Cellenion SASU, Lyon, France; ²Scienion GmbH, Berlin, Germany

Overview

Here, we present an automated cellular aggregate sorter and dispenser (spheroONE[®]) that enables rapid and standardized 3D spheroid, organoid and tumoroid sample preparation.



Introduction

Complex three-dimensional (3D) in vitro models, in particular spheroids, tumoroids and organoids, offer unprecedented means to yield highly predictable models of healthy and diseased tissues and organs.

- Cancer spheroid mimic in vitro several hallmarks of solid tumor, i.e.:
- Drua resistance
- Differential gene expression
- Hypoxia/necrotic core
- Stem cell-derived organoids replicate native tissue/organ
- Phenotype
- Morphology
- Physiological function
- More predictive models for drug screening and basic research
- Reduced use of animal models
- 1. Bulk spheroid formation & culture
- 2. Collect spheroid suspension
- 3. Automated spheroid isolation/sorting
- 4. Incubate 5. Add drugs to
- individual wells 6. Perform your readouts



E.C. Costa et al., Biotech Adv, 2016

Cellenion SASU 60 Avenue Rockefeller Bioserra 2, 69008 LYON France Tel: +33 986 48 70 70 contact@cellenion.com www.cellenion.com

SCIENION GmbH Volmerstr. 7b D-12489 Berlin Germany Tel: +49 (0)30 6392 1700 support@scienion.com www.scienion.com

SCIENION US, Inc 4405 E. Baseline Road Suite #123 Phoenix, AZ. 85042 United States Tel: +1 (888) 988-3842 USsalessupport@scienion.com

Method



Figure 1: Image-based automated single large-particle isolation. a. Picture and schematic representation of the Nano dispensing capillary (NDC) in front of the camera. A Mapping procedure, which tracks objects inside the capillary as droplets are continuously dispensed allows to determine empirically the Ejection Zone (EZ) (ca. Area corresponding to the volume of the next drop). A sedimentation Zone (SZ) is added by the software to account for particle sedimentation. b. Schematic representation of automated imagebased particle isolation, ca. (i) only when a single particle is detected does the robot move on top of the next target well, ensuring a single particle is dispensed in it, and (ii) if no particle or multiple particles are detected, the next drop is dispensed directly into a recovery vial lying directly below the NDC.

Cell culture & bulk spheroids formation

HEK (or HEK-GFP) cells were cultured in DMEM media supplemented with 10% FBS and antibiotics (Penicillin/streptavidin). HEK spheroid were prepared by liquid overlay cell culture on non-adherent surface under constant agitation for 3-7 days. Culture media was exchanged every second day. Before processing on the spheroONE, spheroid suspension. Pre-differentiated intestine organoids (day 4) were prepared on Gri3D plate and kindly provided by SunBiosciences.

SCIENION (UK) Ltd 2000, Lakeside North Harbour Western Road, Portsmouth PO6 3EN United Kingdom +44 (0)7483 388 271 +44 (0)23 9323 3603 support@scienion.com

Automated spheroid isolation and dispense



Gentle dispensing maintains aggregate integrity and viability

Cellular aggregate sorting by size



Advanced sorting



Sort aggregates based on diffraction pattern to specifically select cysts and lumenized organoids (Intestine) ▶ Re-encapsulate aggregates in ECM using pre-filled target plate placed on chilled (4°C) plate holder

Conclusion

Accuracy

Up to 100% single spheroid per well

Homogeneity

• User-defined sorting of spheroids by size and shape

each wells

Advanced sortina

- on diffraction pattern
- Compatibility with organoid workflow

Versatility

- Open-platform
- Any standard well plates (i.e. 96, 384) or specialty/custom designed labware
- Temperature of target labware and environmental humidity control



