

MANUAL

proteoCHIP LF 48

for the Label Free
Single Cell Proteomics workflow
using the cellenONE®




Table of Contents

1. Introduction.....	3
2. Materials	4
3. cellenONE setup	5
3.1. Nozzle setup and Drop parameters setup	5
3.2. Humidity and temperature setup	6
3.3. Target setup.....	8
3.4. Field setup	9
3.5. Drop Setup & Start Dispensing.....	10
3.6. proteoCHIPs inside the cellenONE	12
4. Step 1: Master Mix dispensing.....	13
4.1. Probe, Run and Target selection	13
4.2. Field selection	14
4.3. Start the Run.....	16
5. Step 2: Cell isolation	17
5.1. Probe, Run and Target selection	17
5.2. Field Selection	17
5.3. cellenONE module	19
5.3.1. Setting the detection parameters.....	21
5.3.2. Setting the ejection boundary	23
5.3.3. Setting the isolation parameters	24
5.3.4. Checking the setup of your experiment.....	25
6. Step 3: Incubation at 50 °C	25
6.1. Probe, Run and Target selection	26
6.2. Field Selection	26
6.3 Start the Run	28
7. Step 4: Cooldown	29
7.1. Probe, run and target selection.....	29
7.2. Start the run	30
8. Step 5: Dilution.....	31
8.1. Manual dilution.....	31
8.2 cellenONE automated dilution	31
8.2.1. Probe, Run and Target selection.....	31
8.2.2. Field Selection	32
8.2.3 Start the run	33
9. Sample transfer	34

1. Introduction

The cellenONE software monitors cells inside the piezo dispensing capillary before dispensing. All steps associated with single cell isolation and dispensing experiments are presented in the Software User Guide.

Note: It is highly recommended to read the Software User Guide before performing any experiment with the cellenONE.

This Manual will guide you through all steps of the workflow dedicated to single cell proteomics label-free experiments using proteoCHIP LF 48.



IBSCI™: Image-Based Single Cell Isolation

IBSCI™ is a high resolution-based technology consisting of automated image acquisition, processing and advanced algorithms to automatically isolate single cells from a cell suspension.

FIBSCI™: Fluorescence Image-Based Single Cell Isolation

FIBSCI™ utilizes multichannel fluorescence-based images to identify and isolate single cells of interest from subpopulations in mixed cell samples based on the presence, absence and intensity of a fluorescent signal.

2. Materials

Part number	Description	
P-20-CL	cellenONE PDC L Piezo Dispensing Capillary (Fixed drop volume between 450-600 pl) 4x manifold	
P-20-CM	cellenONE PDC M Piezo Dispensing Capillary (Fixed drop volume between 300-450 pl) 4x manifold	
C-PLF-48-3	proteoCHIP_LF_48 The proteoCHIP LF 48 Set containing: <ul style="list-style-type: none"> • 3 (three) proteoCHIP LF 48 	
	The proteoCHIP LF 48 is a pyramidal-shaped bottom nanowell chip, 48 nanowells, prefilled with oil, enabling label-free single cell proteomics sample preparation.	
	proteoCHIP_LF_48 cellenONE Holder Holder to accommodate up to 3 proteoCHIP_LF_48 on the deck of the cellenONE throughout sample preparation (cell isolation, reagent dispensing and incubation steps). <i>(1 piece/pack)</i>	
C-PLF-48-CHT	proteoCHIP_LF_48 cellenONE Holder Top Top to maintain the proteoCHIP_LF_48 inside the cellenONE Holder on the deck of the cellenONE throughout sample preparation (cell isolation, reagent dispensing and incubation steps). <i>(1 piece/pack)</i>	

3. cellenONE setup

3.1. Nozzle setup and Drop parameters setup

Before starting your experiment, you will need to make sure that your PDC is correctly positioned and produces stable droplets.

At the beginning of the priming, the PDC should be centred with the red cross. You can align your PDC using:

- Nozzle Setup >> Nozzle Offset >> Up, Down, Back, Front, Left, Right controls (Figure 1)

After the priming enter the Voltage and Pulse parameters (indicated on the PDC box), then click on:

- Set Nozzle Parameters
- Start Continuous Dispensing

If needed, adjust the Voltage and Pulse until you have a stable drop and then press:

- Set Nozzle Parameters



Figure 1. Drop optimisation

3.2. Humidity and temperature setup

Throughout the workflow, the humidity and temperature need to be set to a specific value to avoid evaporation and minimize condensation.

Note. The runs already contain specific tasks controlling humidity and temperature for each run, but you should verify that the dew point temperature (DP) is above the temperature in the enclosure (BT) to avoid evaporation ($BT < DP$). We recommend setting -1 °C dew point correction (see Figure 4).

The following describes how to change humidity/temperature “manually” if you wish to change it.

- Humidity setup

Main >> File >> Humidity Control >> Popup window >> Set humidity >> Click on “Set” (Figure 2).




Figure 2. Humidity control

- Temperature setup

Main >> File >> Cooling Unit Control >> Popup window >> Set the temperature to a fixed control or at Dew point and turn the cooling unit ON.

Please also make sure that the Dew Point Correction Value (Cooling Unit Control >> Popup window >> Settings tab) is set to -1 (Figure 4).




Figure 3. Temperature control




Figure 4: Dew Point Correction Value

Note. During a long cell isolation run, the -1 °C Dew Point Correction Value can result in some condensation on the plate. However, this is not affecting the Master Mix efficiency neither is it leading to cross contamination.

Note. When doing a label free experiment, it is best to start by setting the temperature to 10 °C.

3.3. Target setup

Note. To teach the proteoCHIP_LF_48, please refer to the “proteoCHIP_User Manual_proteoCHIP_LF_48_Teaching” manual.

To display your target, go to:

- Target Setup >> Target

For the proteoCHIP LF, you will always have the following configuration in Target Setup:

- Start point: 0, 0
- Field Gap: 0,0
- Number of Fields: X=1, Y=1
- Spot per field: X=12, Y=4
- Spot pitch: X=4500 μm , Y=4500 μm




Figure 5. Target setup

3.4. Field setup

The “Field Setup” sub-tab is used to define:

- From which well the sample will be aspirated.
→ *ONLY if the run contains a “TakeProbe” task (This is the case in the LF protocol).*
- Number of drops/volume to be dispensed in which “spot” of the Target.

Note. To change from Volume to drops and vice versa:

- Target Setup >> Field Setup >> “Field” at the top >> Volume entry

Note. For every run, fields are saved in the software.

You can edit these Fields using “Edit Field Table” to change:

- From which well the sample will be aspirated.
- Number of drops/volume to be dispensed in which “spot” of the Target.




Figure 6. Field setup

Note. To assure that the exact suggested volume of each reagent is dispensed we recommend performing a Drop volume to calculate the corresponding number of drops before each liquid dispensing run.

This should be done before dispensing the Master Mix, H₂O in the rehydration steps and 0.1% TFA/FA.

To modify the number of drops/volume in the Field setup refer to Figure 7 and change the number of drops/volume by clicking on the Drops subtab.




Figure 7. Edit Field Table

3.5. Drop Setup & Start Dispensing

Each dispensing run is built with the same logic, once the run starts (Figure 8):

Several pop-up messages will prompt you to load the correct volume of the desired solution in the well specified in the Target Setup section.

- The sample should be loaded in a chosen well of the 384 source plate (The well is chosen in the Field Setup tab, see paragraph 3.4).
- The volume to be aspirated is predefined inside the Run by the “TakeProbe” task. (Required volume will be mentioned in the message that will pop up at the beginning of the run).

- The capillary (PDC) moves to the chosen well, aspirates the chosen volume and moves back to do a wash cycle and take a picture of the drop.
- The software will automatically check that the dispensing parameters are correct by performing a “DropCheck”. Should the system fail to detect a suitable drop, the run will pause and the user will be able to modify the parameters to restore a stable droplet generation.
- The user must check that the dispensing parameters (Voltage and Pulse) (Figure 8 (1)) are correct by pressing the “Continuous Dispensing” button (Figure 8 (2)) until a drop is stable without any satellite (if any).
- The user must save the new settings if the Voltage and Pulse were modified (Figure 8 (3)) and press the “Drop Check” button (Figure 8 (4)) and check if the drop is identified by the software (circled in red) (Figure 8 (5)).
- If needed, you can also perform troubleshooting tasks (AirEx, Flush_PDC, etc.) and uptake more sample.
- Once the droplet generation is stabilized, click “Continue” to resume the run (Figure 8 (1)).




Figure 8. Drop Setup and start dispensing

3.6. proteoCHIPs inside the cellenONE




Figure 9. proteoCHIP LF inside the cellenONE holder

- Place the proteoCHIP(s)_LF_48 inside the cellenONE holder (ref. C-PLF-48-CHB) (Figure 9).

Note. If processing only one chip, insert it on the spot closest to you.

- Place the cellenONE holder top (ref. C-PLF-48-CHT) on top of the holder and insert the assembled unit on the target holder inside the instrument.

Note. The experiment is a standardized workflow for only one proteoCHIP LF 48. If you are doing more than one proteoCHIP LF 48 in parallel, you will need to repeat run 5 for each proteoCHIP LF 48.

e.g., If you are running 3 proteoCHIPS you will need to do:

- 1X the run 1
- 1X the run 2
- 1X the run 3
- 1X the run 4
- 3X the run 5

4. Step 1: Master Mix dispensing

Note.

We recommend using this Master Mix (Lysis + Digestion):

0.2% DDM, 10 ng/ μ L Enzyme(s), 100 mM TEAB

For each experiment 100 μ L of Master Mix is needed.

To prepare 100 μ L of Master Mix add the reagents in the following order:

10 μ L 1M TEAB + 60 μ L H₂O + 20 μ L 1% DDM + 10 μ L 100 ng/ μ L Enzyme(s)

Do not vortex the mixture, this will introduce air bubbles. Mix by gently pipetting the solution up and down.

Tip: Aliquot 20 μ L 1% DDM and 10 μ L 100 ng/ μ L Enzyme(s) at -20 °C for short term usage.

For long term storage of enzymes, we recommend storage at -80 °C.

TEAB should be stored at 4 °C.

Note.

We recommend using MS-grade enzymes (e.g., Trypsin Gold, Promega or Trypsin/LysC Mix, Thermo). For enzyme mixtures we recommend to also use a 10 ng/ μ L mixture of the enzymes.

Note.

For the Master Mix, you will probably need higher Voltage (+5-10%) to achieve an optimal drop formation. You can also optimize that during the run by clicking Nozzle Setup.

4.1. Probe, Run and Target selection

The Master Mix dispensing happens at 8 (+/- 2) °C and at humidity set to 45%.

In the Main tab of the software, select the probe, run and target to be used (Figure 10).

- Probe: "MTP384"
- Run "LF_48_Run_1_MasterMix"
- Target: "proteoCHIP_LF_48"

This run contains an intermediate task that checks the stability of the drop after the dispensing of each row of the proteoCHIP (Figure 15). If the drop is not stable, you will have the opportunity to optimize the parameters, as explained in paragraph 3.5.




Figure 10. Run 1 Master mix dispensing

4.2. Field selection

Load the field “LF_48_Run_1_MasterMix” which will spot 400 drops twice (~300 nL) per well.

Note.

As the run will dispense the Master Mix in two rounds, to have the correct end volume (~300 nL), you need to modify the field in volume entry to dispense 150 nl.

- To load the field:

“Target Setup” >> “Target” >> “Load Field” >> Select the field corresponding to the run (Figure 11) and you will have the field corresponding to Figure 12.




Figure 11. Load Field for MasterMix dispensing




Figure 12. Field Run 1

4.3. Start the Run

Note. Load 60 μL of the Master Mix in the selected well.

Once the proteoCHiPs are loaded in the cellenONE as seen in Figure 9, check that you selected the right probe, run and target in the run window (Figure 13(1)). Start the run (Figure 13(2)), name the experiment (Figure 13(3)) and click OK (Figure 13(4)). Select the number of proteoCHiPs that you want to dispense (Figure 13(5)) and click “OK” (Figure 13(6)).

A few pop-up messages will guide you before starting to spot (Figure 14).




Figure 13. Starting the MasterMix dispensing run




Figure 14. Pop-up messages

Note. ALWAYS use at least the target 1 so the “AutoDropDetection” and “OpenNozzleSetup” works properly (first proteoCHIP as shown above).




Figure 15. Autodrop check

5. Step 2: Cell isolation

Prepare your cells at a concentration of 100-200 cells/ μ L in degassed PBS.

5.1. Probe, Run and Target selection

In the Main tab of the software, select the probe, run and target to be used (Figure 16).

- Probe: “MTP384”
- Run “LF_48_Run_1_MasterMix”
- Target: “proteoCHIP_LF_48”




Figure 16. Run 2 Single cell isolation

5.2. Field Selection

Click on the 3rd tab, “Target Setup”, and the 1st sub-tab, “Target”. Target is the same for the whole workflow, but the field setup will need to be changed at every step.

Note. To isolate and dispense single cells, positions defined in the field setup must contain only one drop per well.

Load the field “LF_48_Run_2_Cells” which will spot 1 drop (=1 cell) per well. To load the field, click in “Target” (Figure 17) >> “Load Field” (Figure 17), select the field corresponding to the run (Figure 17) and you will have the field corresponding to the Figure 18.




Figure 17. Load Field for cell isolation




Figure 18. Field Setup for cell isolation

Note. If you wish to dispense more than 1 single cell, the following syntax must be used: 1, 1 (for isolation of two single cells); 1, 1, 1 (for isolation of three single cells); 1,1,1,1 (for isolation of four single cells), etc. You can do so by clicking multiple times (n) on the well where you want to have 1+n cells isolated.

5.3. cellenONE module

Note. Load > 30 μL of cell suspension (100-200 cells/ μL) into a chosen well of the 384-well probe plate.

For this run you will need to **open the cellenONE window** and start the run in this same window. Once the **cellenONE module** is opened, it **should not be closed** before the end of the run (nevertheless the window can be minimized).

Please refer to the Software Reference Guide Chapter 7 “Single cell isolation” to understand all the options and features for single cell dispensing.

Select the “Nozzle Setup” tab (Figure 19 (1)) and press on the cellenONE logo (Figure 19) to open the cellenONE module for single cell isolation.

The cell isolation run is still performed at 10 °C and ~45% humidity.




Figure 19. cellenONE module for single cell isolation (right window)

Before starting, define the folder location and name your experiment by pressing the following

button 

All files associated with your experiment will be saved in this location. By default, the data is saved in the last folder chosen and named `samplename_date`.

Start the run by clicking on the .



Please note that you can always pause the run with the Nozzle setup shortcut . You will be able to modify the parameters if wanted.

Select the proteoCHiPs that you want to dispense in (Figure 20) and click “OK” (Figure 20(4)).




Figure 20. Starting the run for single cell isolation




Figure 21. Pop-up messages

After uptaking the desired volume of the cell sample, the PDC comes back to the camera station so you can define the cell parameters used for cell isolation.

5.3.1. Setting the detection parameters

The “Main” tab of cellenONE module (Figure 22) displays images of the region of interest of the PDC in which cells are detected. During dispensing, cells are travelling from right (upper part of the nozzle) to left (nozzle tip). Information below the image (i.e., Diameter, Elongation, Circularity, and Intensity) is given for the leftmost detected cell, i.e., the closest to the nozzle tip.




Figure 22. Main tab of cellenONE® module.

Hint. The elongation factor can be used to effectively differentiate two aggregated cells from a single cell as these will have different elongation factors.

Round single cells: <1.6; Doublets >2.5

This section is used to define detection and isolation parameters (Figure 23). These parameters are crucial to successfully isolate and dispense single cells.




Figure 23. Detection and isolation parameters

Detection parameters: These parameters will determine what the software will detect and consider as an object (e.g., cell or particle). As a result, these parameters should be “large” enough to ensure that all cells and cell clusters are detected as objects.

- **Min diameter:** detected diameter of the smallest object that should be detected. An object with a diameter lower than the entered value will **not be detected**.
- **Max diameter:** detected diameter of the largest object that should be detected. Beware that cell aggregates will typically be detected as a single large object. An object with a diameter bigger than the entered value will **not be detected**.
- **Elongation:** an object with an elongation above (for MAX) or below (for MIN) the entered value will **not be detected**.

HINT 1 The detected diameter is typically about 10% larger than real diameter.

HINT 2 A detected object is surrounded by a green or a red circle:

- a **green circle** means that the object fulfils the isolation parameters and would be isolated.
- a **red circle** means that the object does not fulfil the isolation parameters and would not be isolated.

HINT 3 To start, set wide detection parameters and then, restrict them if needed.

- Diameter min: 5 µm
- Diameter max: 100 µm
- MAX Elongation: 4

After loading the sample, do a background:



Get background button: Takes a reference background image. This button should only be pressed when the nozzle is free of cells and/or particles and when the PDC is in front of the camera.

Make several drops:



Test droplet button: Ejects a single droplet.

Check visually by eye that every particle inside the PDC is detected (i.e., that the particle is circled by a red or green circle).

5.3.2. Setting the ejection boundary

Once detection parameters are set, do a mapping:



Start mapping button: Starts the mapping process.

This process maps the position of cells inside the PDC during the dispensing of droplets in front of the camera. This will determine **the size of the ejection zone** (Zone 1 on Figure 22). Pressing this button during mapping will interrupt the mapping process.

At the end of the mapping process, the ejection boundary is set automatically according to the blue and green dot distribution inside the PDC (Figure 24). Cells that were detected in the location of **blue dots** were subsequently detected again inside the PDC upon next droplet generation. Whereas **green dots** correspond to the last seen location of a cells within the PDC (i.e., upon next droplet generation these cells were ejected from the PDC). Once done, click “Save & close”.




Figure 24. Mapping

Hint. At the end of the mapping process, the ejection boundary is set automatically according to the distribution of the blue and green dots inside the PDC. To **adjust the ejection zone manually, move the slider** by clicking and dragging the positioner on the ruler below the image and **press “Set” to validate.**

5.3.3. Setting the isolation parameters

Once Mapping is done, **set up your isolation parameters.**

Isolation parameters: These parameters can be determined from the characteristics of the cells shown in the mapping report (Figure 24) and adjusted as needed to include only cells of interest. These values are critical to ensure satisfactory single cell isolation. All objects that fulfil the isolation parameters and single cell condition will be isolated.

- **Min diameter:** detected diameter of the smallest object that should be isolated.
- **Max diameter:** detected diameter of the largest object that should be isolated.
- **Elongation:** an object with an elongation above (for MAX) or below (for MIN) the entered value will **not be isolated.**

Three options are available here:

- Draw a gate right after the mapping from the scatter plot that was just performed or enter values manually into the isolation parameters based on the mapping.
- Open the mapping (in the folder you chose earlier, folder name.mapping) and interpret the different scatter plots.
- Perform an analysis of your sample.



Analyze sample button: Analyzes 100 cells from your sample.

This can be used to estimate the heterogeneity of your sample. At the end of the analysis, a window will pop up with scatter plots (Figure 25). The number of analysed cells can be set in the “Advanced Settings” tab. You can do a gate to choose new values for isolation parameters.




Figure 25. Scatter Plot after Analysis

5.3.4. Checking the setup of your experiment

Check the ejection boundary:



Test droplet button: Press it manually multiple times to verify the ejection boundary. When a particle is inside the ejection zone, it should be ejected within the next droplet when pressing the droplet button once more.

Check the isolation parameters:

cells (or particles) of interest should be circled with a **green circle** and the other particles with a **red circle**.

Note. Continuous drop generation can be activated by simultaneously clicking on the left mouse button and Shift button and is interrupted by clicking once again on the test droplet button.

Once you are satisfied with your parameters you can click Continue to start dispensing your cells (Figure 26).




Figure 26. Starting cell isolation Run

6. Step 3: Incubation at 50 °C

Single cells have been isolated in wells containing lysis and digestion buffer. The next step is to perform the incubation to:

- Lyse the cells,
- Digest the protein into peptide fragments.

This step can be performed directly on the deck of the cellenONE. This incubation step will last 1h30, at a temperature of 50 °C. In addition to an increased ambient humidity, the PDC will regularly dispense deionized water in each well to maintain a constant water content and create a mechanical mixing phenomenon to help achieving homogeneous cell lysis and protein digestion.

6.1. Probe, Run and Target selection

In the Main tab of the software, select the probe, run and target to be used (Figure 27).

- Probe: “MTP384”
- Run “LF_48_Run_3_Incubation_50C”
- Target: “proteoCHIP_LF_48”




Figure 27. Run 3 Incubation

6.2. Field Selection

Click on the 3rd tab, “Target Setup”, and the 1st sub-tab, “Target”. Target is the same for the whole workflow, but the field setup will need to be changed at every step.

Load the field “LF_48_Run_3&4_Incubation&Cooldown” which will spot 315 drops (~120 nL) per well. To load the field, click in “Target” Figure 17 >> “Load Field”, select the field corresponding to the run (Figure 28) and you will have the field corresponding to the Figure 29.




Figure 28. Load field for incubation




Figure 29. Field setup for incubation

6.3 Start the Run

Start the run and **select all 3 chips**, then click OK (Figure 30).

Warning. Please make sure that the 3 proteoCHiPs are selected regardless of the number of chips processed.




Figure 30. Start the run for Incubation

Once the Run has started, humidity and temperature are automatically set. A pop-up message will be displayed, click “OK” to start the 1h30 incubation (Figure 31).




Figure 31. Incubation messages

After 1h30, once you are ready to proceed to the next step, click the red “Stop” button. Dismiss the pop-ups by clicking “Exit Run” and then “Partially spotted” (Figure 32).



Figure 32. Exit the run

7. Step 4: Cooldown

Following incubation at 50 °C, we need to bring the plate temperature back to 20 °C. During this process, the system will keep on rehydrating the wells continuously. The cooldown step takes approximately 30 minutes to complete.

This run needs to be launched **immediately** after stopping the previous run, as in the meantime rehydration cycles are interrupted, but the temperature setting is still 50 °C. Not starting the Cooldown run immediately can result in sample evaporation.

7.1. Probe, run and target selection

In the Main tab of the software, select the probe, run and target to be used (Figure 33).

- Probe: “MTP384”
- Run “LF_48_Run_4_Cooldown_20C”
- Target: “proteoCHIP_LF_48”

No changes need to be made in the field setup section, so you can proceed to start the run immediately.




Figure 33. Run 4 Cooldown to 20 °C

7.2. Start the run

Once the run has been selected, click “Start Run” in the Run tab. Select all 3 chips regardless of the number of loaded chips and click “OK”.

The run will finish once the temperature of 20 °C is reached. This should take approximately 30 minutes.




Figure 34. Starting the run for Cooldown




Figure 35. Pop-up messages




Figure 36. Incubation completed

8. Step 5: Dilution

The dilution step is used to adjust the end volume of the final sample that is going to be injected in the LC. The final volume should be around 3.5 μL to be correctly injected into the HPLC system loop before the MS analysis. To achieve that, 0.1% FA or 0.1% TFA is added to each well of the chip.

This step can either be performed manually or by an automated cellenONE run.

8.1. Manual dilution

Please make sure that the oil is liquid (temperature should be above 18 $^{\circ}\text{C}$). Pipette 3.2 μL of 0.1% FA or TFA in each well of the chip. Optional: Use a multichannel pipette.




Figure 37. Manual dilution using a multichannel pipette

8.2 cellenONE automated dilution

8.2.1. Probe, Run and Target selection

In the Main tab of the software, select the probe, run and target to be used (Figure 38).

- Probe: “cellenWASH” / “cellenWASH_TALL”
- Run: “LF_48_Run_5_Dilution_1chip”
- Target: “proteoCHIP_LF_48”




Figure 40. Target setup Dilution

8.2.3 Start the run

Once you chose the right Probe, Run, Target, and you loaded the dedicated field, click “Start Run” in the Run tab (Figure 41). A pop-up window will prompt you to insert a cellenVIAL or a pre-taught PCR tube (cellenWASH_TALL) loaded with 250 µL of 0.1% TFA or FA in the cellenWASH1 position (Figure 41). Once this is done, click “OK” to continue. If the cellenWASH_TALL is not taught on your machine and you would like to use a PCR tube, please contact your Field Application Specialist.





Figure 41. Starting the Dilution run

9. Sample transfer

Place the holder containing the chip(s) on ice so the temperature decreases, and the oil becomes solid (temperature below 18 °C) (Figure 42).

Once the oil is solid (after ~5 minutes), it is easy to pipette out 3.5 µL of the sample droplet. The solidified oil will keep sticking onto the chip and will not be transferred with the sample to be analysed.

You can now transfer your samples to the desired injection plate/vials, etc.



Figure 42. Solid oil in proteoCHIP LF 48

Note. For the moment this step is performed manually. We are currently working on implementing direct injection from the chip to most of the common autosamplers available. As soon as this feature is available, we will let you know!