

Drug response profiling of individual primary colorectal cancer tumoroids using a novel automation workflow and AI-assisted image analysis



Anthony Eritano¹, Dylan Bryda², Andreea Merloiu³, Gerri Ortiz¹, Ekaterina Nikolov³, Barbara McIntosh¹, Seyoum Ayeahunie², and Evan F Cromwell³

¹Scienion, Inc. 1231 W. Warner Rd, Suite 103, Tempe AZ 85284; ²MatTek Inc, 200 Homer Ave, Ashland MA 01721; ³PFI, Inc. 875 Cowan Rd, Suite B, Burlingame CA 94010

INTRODUCTION

The use of human-relevant three-dimensional *in vitro* cell models have gained popularity as they better recapitulate key aspects of human tissues and tumor microenvironments. Generating organoids, spheroids, and tumoroids have progressed to where scientists can recapitulate most human organs and cancer types *ex vivo*.^{1,2} These models are being used in many aspects of drug discovery and development, but adoption has been limited due in part to challenges related to sample handling, assay development and the need for optimized instrumentation.³ Here we present a novel automated workflow that addresses those issues. We have combined the use of spheroONE[®] for automated sorting, isolation and dispensing of tumoroids and Pu-MA System[®] EC with low attachment flowchips for performing viability assays and biomarker detection.



The spheroONE, a large-particle sorter and dispenser, was optimized to dispense size-selected tumoroids into the protected sample chamber of Pu-MA System flowchips.^{4,5} The Pu-MA System EC performed compound additions and staining using automated fluid transfers. Drug treated and untreated tumoroids were incubated for 48 hr in flowchips in a standard incubator to maximize efficiency of Pu-MA System. The samples were stained in the flowchips with viability dyes, imaged with ECHO microscope, and analyzed with an AI-assisted program. We used this workflow to analyze proliferation and viability of HCT116 colorectal cancer (CRC) spheroids and patient-derived CRC tumoroids. This platform is a valuable tool in a wide range of research areas including disease modeling, drug discovery and personalized medicine.

AUTOMATED REAGENT EXCHANGE

The Pu-MA System EC and 3D Flowchip features:

- Automated media exchanges occur with cells in protected chamber
- Supernatants can be collected to monitor cell secretion
- Spheroids can be stained & imaged in the flowchips
- Assay protocols can be edited via the Pu-MA System Software




Figure 1. Pu-MA System EC workflow allowing media exchange, drug treatment, and immunostaining of individual spheroids. The Pu-MA System EC system workflow allows for the processing of 32 individual spheroids inside an environmentally controlled tabletop device. Precise microfluidic channels allow for the addition, removal, and recovery of up to seven different reagents/samples. Spheroids are contained within a microcavity to prevent shear damage during fluid exchange. High clarity bottom allows for excellent imaging with an inverted scope.

SPHEROID ISOLATION AND DISPENSE

spheroONE features:

- ≥95% efficiency for automated dispensing of patient derived CRC tumoroids using morphology (size and shape)
- Size based sorting from 100 to 600 μm diameter
- Gentle dispensing technology maintains the integrity and viability of fragile cellular aggregates.
- Direct Visual Inspection of the sample along with full image record




Figure 2. SpheroONE workflow to isolate spheroids. (A) Cells are grown and placed into spheroid forming microcavity plate. (B) Spheroids are collected and transferred to a sample reservoir. (C) Sample and desired plate are loaded into the spheroONE. (D) Using Image-Based Isolation Software, the spheroONE images a dispensing capillary. The capillary is divided into three zones: The bubble Zone (orange) represents the volume displaced during dispense. The Ejection Zone (green) represents volume of media that will be dispensed in next drop. The Sedimentation Zone (blue) represents the region where spheroids settle and can potentially be dispensed.

PATIENT DERIVED CRC TUMOROID WORKFLOW




Figure 3. Patient derived primary tumors were obtained from Next Oncology (San Antonio, TX) as primary or secondary passages. Further passaging and expansion of the tumoroids were performed at MatTek. Cancer epithelial cells were expanded from frozen stock and then seeded into multicavity plates (Heidolph Sphericalplate 5D 24-well plate) with a seeding density of approximately 600 cells per microcavity. Tumoroids were allowed to form over a 10 - 13 day period before harvesting and loading into the spheroONE for automated dispense into Pu-MA System flowchips.

VIABILITY ASSAY SETUP & WORKFLOW

Assay reagents are loaded into flowchip wells and then fluid exchanges occur with the sample chamber according to pre-determined assay protocols. Tumoroids were protected during these operations and approximately 95% of fluid is exchanged. This offers efficient washing and minimizes compound carry-over. The assay protocol here included compound incubation, viability staining, and wash steps. Tumoroids were characterized by automated high resolution transmitted light and fluorescence imaging.




Figure 4. Workflow for viability assays using P-D CRC tumoroids. (A,B) Multiple plates were loaded with tumoroids and compounds were dispensed in serial dilutions (e.g., 32 samples per plate 3-fold serial dilution, 8 conc., n = 4). Flowchips were placed in incubators using holders with water reservoirs to limit evaporation. After 48 hr incubation, staining solution (CyQuant Green + EthD-1 or PI, Invitrogen) was loaded and transferred into the sample wells using Pu-MA System. Imaging was done with either ECHO Revolution (10X, widefield, ECHO) or CQ1 (10X, confocal, Yokogawa).

AI-ASSISTED VIABILITY ANALYSIS




Figure 5. Workflow of AI-assisted analysis of fluorescently stained spheroids for viability assessment. We have developed a comprehensive method for automated spheroid detection and analysis of fluorescent images, employing state-of-the-art deep learning techniques. The approach combines the power of transfer learning using a pre-trained VGG16 regional convolutional neural network (R-CNN) with custom-built fully connected layers for precise spheroid detection.

- Images undergo preprocessing operations to ensure compatibility with the VGG16 model (standardization of image size and intensity normalization).
- Pre-trained VGG16 R-CNN is employed to extract features from the images.
- Custom-built fully connected layers learn discriminative features through iterative training specific to spheroid detection.
- The model, facilitated by the Adam optimizer and binary cross-entropy loss function, measures spheroid parameters and produces data outputs.

Parameter	Value
Total Area	3449 px
Total Mean Intensity	68549
Total Intensity per Area	872101

References

- Matossian, M.D. et al (2019) BMC Cancer. 19(1): 205
- Matossian, M.D. et al (2021) Clin. & Trans. Oncology 24: 127-144
- Sirenko, O. et al (2015) ADDT 13, 402
- Cromwell, E.F. et al (2021) SLAS Tech. 26(3): 237
- Cromwell, E.F. et al (2022) SLAS Disc. 27: 191.
- Khaleel, S.A. et al (2016) Sci Reports 6, 36855
- Narvi, E. et al, (2018) Sci Reports 8, 16579
- Richard, S.M. et al, (2015) 11, 336

AUTOMATED FLOWCHIP LOADING

An important design criteria of the spheroONE is to not adversely affect spheroid viability. To characterize this, sister spheroids were dispensed into flowchips via manual pipetting or spheroONE. Sample wells were filled with 20 μL of media and incubated overnight. Media was replaced with viability staining solution using Pu-MA System and subsequently imaged, and then spheroids were imaged for Live and Dead cells. Equivalent phenotypes and viability were observed between two methods with % Live cells > 90%.




Figure 6. (A) Comparison of spheroid viability using manual and spheroONE dispense in Pu-MA System flowchips. (B) Images captured by spheroONE showing isolation of tumoroids and position in flowchips. (C) Analysis of LIVE/DEAD staining showing survivability of manual and spheroONE isolated spheroids

	Percent Live	
	Average*	StDev
Manual	94%	3.8%
spheroONE	97%	1.3%

*n=5

- Imaged based isolation ensures only spheroids are dispensed
- Post-dispense image shows that correct spatial dispense
- Visualization under microscope confirms isolated of spheroid into center microcavity

HCT116 COMPOUND RESPONSE

Spheroids were treated for 48 hours with compounds doxorubicin and oxaliplatin then stained with viability dye. Spheroids were imaged in the flowchips and analyzed for number of Live and Dead cells. Spheroid viability results were found to be consistent with published response of the compounds for 2D cell cultures (0.96 μM).⁶





Figure 7. (A) Composite MIP images of HCT116 spheroids treated with different concentrations of doxorubicin and oxaliplatin. Viability-response curve for incubation with compound. Image analysis was done with AI-assisted process. (B) Live (Green) and Dead (Red) channels were analyzed independently. Live Cells value was determined from normalized integrated intensities: Live = IntGr/(IntGr + IntRed) Data was fit with 4P function. (Error bars = +/- 1 SD, n=4)


	Doxorubicin	Oxaliplatin
IC ₅₀	0.75	11.8
Slope	2.9	2.44

- Doxorubicin⁶: 0.96 μM
- Oxaliplatin⁷: ~5 μM

HCT116 VS PATIENT-DERIVED CRC TUMOROID



- B Chemotherapeutics**
- Doxorubicin**
 - Anthracycline, blocks enzyme topoisomerase 2. Widely used chemo agent susceptible to drug resistance.⁶
 - P-D CRC and HCT116 show similar sensitivity
 - Oxaliplatin**
 - Ligand-Pt complex causes DNA damage, induces apoptosis.^{7,8}
 - P-D CRC show less sensitivity than HCT116
 - 5-Fluorouracil**
 - Precursor of dTTP and UTP. Interferes with both DNA and RNA metabolism affecting DNA repair and DNA or RNA synthesis.⁸
 - P-D CRC ~50x more sensitive than HCT116
 - Cetuximab**
 - EGFR mAb for RAS wild-type metastatic CRC.⁷
 - P-D CRC sensitivity seen at 68 ug/mL.
 - HCT116 not sensitive



Compound	IC ₅₀ (μM)	Slope
Doxorubicin	0.49	1.3
Oxaliplatin	55	0.57
5-Fluorouracil	1130	6.2
Cetuximab	n.r.	

Figure 8. (A) Representative images of P-D CRC tumoroids and HCT116 spheroids after incubation with compounds for 48 hrs. Spheroids were stained with CyQuant Green (Live) and EthD-1 or PI (Dead). (B) List of chemotherapeutic drugs used for initial trials. (C) Dosage curve showing IC₅₀ for the three compounds. Cetuximab showed no response in HCT116 Cells. Imaging was done with either ECHO Revolution (10X, widefield, ECHO) or CQ1 (10X, confocal, Yokogawa).

EGFR INTERNALIZATION DISRUPTED BY CETUXIMAB




Figure 9. HCT116 cells treated with human recombinant EGF and/or Cetuximab (250 ug/mL) overnight then treated with EGF (10 nM) for 20 min. Cells were fixed and stained for EGFR and F-actin. Images were acquired using a CQ1 confocal imaging system (40X Obj.). (B) EGFR internalization was analyzed using CellPathfinder (Yokogawa) and showed a shift from membrane localization (bottom mask) to cytoplasmic localization (top mask) of EGFR upon EGF activation. Cetuximab was observed to interfere with this internalization.

CONCLUSIONS

- We have demonstrated capabilities of a novel automated spheroid and tumoroid assay system that performs complex protocols with 3D cell models.
- Tumoroids are automatically sorted, isolated and dispensed into flowchips for downstream assays providing precise control over size and number.
- Assays and fluid exchanges are performed in a novel microfluidic device that protects the cell models and enhances assay precision and control.
- The ability to analyze spheroids and tumoroids *in situ* in order to capture toxicity information and perform functional assays shows great promise for disease modeling and drug discovery.