


Drug response profiling of individual primary colorectal cancer tumoroids using a novel automation workflow and AI-assisted image analysis

Anthony Eritano^{1*}, Dylan Bryda², Andreea Merloiu³, Alissa Schutte¹, Gerri Ortiz¹, Barbara McIntosh¹, Seyoum Ayeahunie², and Evan F Cromwell³
¹SCIENION, Inc. 1231 W. Warner Rd, Suite 103, Tempe AZ 85284 ²MatTek Inc, 200 Homer Ave, Ashland MA 01721 ³PFI, Inc. 36 Bennett Rd, Redwood City CA 94062
*Presenting author: a.eritano@scienion.com

INTRODUCTION

The use of three-dimensional in vitro cell models have gained popularity as they better recapitulate key aspects of tissues and tumor microenvironments. Generating organoids, spheroids, and tumoroids have progressed to where scientists can emulate most human organs and cancer types *ex vivo*^{1,2}. These models are being used in many aspects of drug discovery and development, but adoption has been limited due in part to challenges related to sample handling, assay standardization and the need for optimized instrumentation³. **Here we present a novel workflow that addresses those issues and analyzed the proliferation and viability of established cell line-derived spheroids (HTC116) and patient-derived CRC (P-D CRC) tumoroids after drug treatment.**



PRECISION LIQUID DISPENSING

- **Accuracy Built-In:** Real-time droplet detection verifies volume and alerts when source liquid runs out.
- **Prevents Cross-Contamination:** Non-contact, pressure-based dispensing eliminates carryover and tip usage.
- **Reagent-Saving Miniaturization:** Dispenses as low as 17.3 nL with <1 µL dead volume to cut reagent costs.
- **Flexible & Compatible:** Supports 17.3 nL–30 µL, handles viscous liquids (up to 43% glycerol), and fits all SBS-format plates, including up to 1536-well plates.
- **Simplified Setup:** Assay Studio software enables protocol creation within minutes - no programming needed

AUTOMATED REAGENT EXCHANGE

- **Automated media exchanges** occur with cells in protected chamber; approximately 95% of fluid is exchanged allowing for efficient washing and minimize compound carryover.
- **Supernatants can be collected** to monitor cell secretion.
- **Spheroids can be stained and imaged** in the flowchips.
- **Assay protocols can be edited** via the Pu-MA System Software.

SINGLE SPHEROID ISOLATION AND DISPENSING

- **Standardization of 3D model size and morphology.**
- **≥95% accuracy** for automated image-based dispensing of patient-derived CRC.
- **Gentle dispensing technology** maintains the integrity and viability of fragile cellular aggregates, from 100 to 600 µm diameter.
- **Direct visual inspection** of the sample along with full image record.

GENERATION OF PATIENT-DERIVED CRC TUMORIDS





Figure 4. Generation of patient-derived CRC tumoroids (P-D CRC). Primary tumors were obtained from Next Oncology (San Antonio, TX) as primary or secondary passages. Further passaging and expansion of the tumoroids were performed at MatTek.

VIABILITY ASSAY DEVELOPMENT



AI-ASSISTED VIABILITY ANALYSIS

- Images undergo **preprocessing** operations to ensure compatibility with the VGG16 model (standardization of image size and intensity normalization).
- **Pre-trained VGG16 regional convolutional neural network (R-CNN)** is employed to extract features from the images.
- **Custom-built fully connected layers** learn discriminative features through iterative training specific to spheroid detection.
- The model, facilitated by the **Adam optimizer** and **binary cross-entropy loss function**, measures spheroid parameters and produces data outputs.

COMPARISON OF MANUAL VS. SPHEROONE SPHEROID ISOLATION






Figure 7. (A) Images captured by spheroONE showing isolation of tumoroids and position in flowchips. **(B)** Comparison of spheroid viability using manual or spheroONE dispensing. Sample wells were filled with 20 µL of media and incubated overnight. Media was replaced with viability staining solution using Pu-MA System and subsequently imaged. **(C)** Comparison of post-isolation viability and single spheroid isolation accuracy.

DOSAGE RESPONSE CURVES FOR IMMORTALIZED CELL LINES



HTC-116 SPHEROIDS vs. PATIENT-DERIVED CRC TUMORIDS



EGFR LOCALIZATION




Figure 10. HTC-116 cells treated with human recombinant EGF and/or Cetuximab. (A) Cells were incubated with Media or Cetuximab (250 µg/mL) overnight then treated with EGF (10 nM) for 20 min. Cells were fixed and stained for EGFR and F-actin. Images were acquired using a CQ1 confocal imaging system (40X Obj). (B) EGF internalization was analyzed using CellPathfinder (Yokogawa)

CONCLUSIONS

- We have demonstrated capabilities of a novel automated 3D cellular model assay system that performs **standardized drug testing**.
- **3D models are automatically standardized and isolated** into flowchips for downstream assays, providing control over size and number.
- Fluid exchanges are performed in a novel microfluidic device that protects the cell models and enhances **assay precision and control**.
- **AI-driven analysis is currently being improved** to recognize and quantify changes in spheroid morphology to allow for further analysis of chemotherapeutic effects on spheroid behavior.
- The ability to analyze spheroids and tumoroids to capture toxicity information and perform functional assays shows great promise for **disease modeling and drug discovery**.

REFERENCES

1. Matossian, M.D. et al (2019) BMC Cancer. 19(1): 205
2. Matossian, M.D. et al (2021) Clin. & Trans. Oncology 24: 127-144
3. Sirenko, O. et al (2015) ADDT 13, 402
4. Cromwell, E.F. et al (2021) SLAS Tech. 26(3): 237
5. Cromwell, E.F. et al (2022) SLAS Disc. 27: 191.
6. Khaleel, S.A. et al (2016) Sci Reports 6, 36855
7. Narvi, E. et al, (2018) Sci Reports 8, 16579
8. Richard, S.M. et al, (2015) 11, 336

SCIENION GmbH, Berlin
Wagner-Régény-Str. 15
D-12489 Berlin

Phone: +49 (0) 30 6392 1700
Fax: +49 (0) 30 6392 1701
Email: support@scienion.com

SCIENION UK Ltd
Building 2000
Lakeside North Harbour
Western Road
Portsmouth
PO6 3EN

Phone: +44 (0) 23 93 233601

SCIENION US Inc.
1231 W. Warner Road, Ste. 103
Tempe, AZ 85284, United States

Phone: +1 (888) 988-3842
USsalesupport@scienion.com

SCIENION Asia, Shanghai
Room 1106B Building 3 No.391
Guiping Road Xuhui District,
Shanghai

Phone: +86 186 2170 9292
support@scienion.com

scienion.com