


CELLENION

A BICO COMPANY

CELLEN ONE HT

True walk-away solution
for high throughput
monoclonal cell line
development




Up to 100% single cell accuracy



True walk-away operation




Assurance of monoclonality



Best clonal outgrowth



Multiple samples to multiple plates



Simply scale-up existing workflows

Key features

- Image-based single cell isolation
- Automated sample loading
- Ultra-gentle acoustic dispensing
- Plate nest shuttle and API
- 4-channel fluorescence imaging


Benefits

- Up to 100% single cell accuracy
- Assurance of monoclonality
- Multiple samples to multiple plates
- True walk-away operation
- Best clonal outgrowth solution on the market
- Simply scale-up existing workflows
- Compatible with 96-, 384-, 1536-wp
- Sort and isolate only cells of interest



Do you want
to see cellenONE
in action?

Book a demo through
our website.



cellenion.com

Contact Us

Cellenion SASU

60 Avenue Rockefeller
Bioserra 2, 69008
LYON France
Tel: +33 986 48 70 70
contact@cellenion.com
www.cellenion.com

SCIENION GmbH

Volmerstr. 7b
D-12489 Berlin
Germany
Tel: +49 (0)30 6392 1700
support@scienion.com
www.scienion.com

SCIENION US, Inc

4405 E. Baseline Road Suite #123
Phoenix, AZ. 85042
United States
Tel: +1 (888) 988-3842
USsalesupport@scienion.com

SCIENION (UK) Ltd

Chichester Enterprise Centre
Terminus Road,
Chichester PO19 8TX
United Kingdom
Tel: +44 (0) 1243 88 71 65
support@scienion.com